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Dispositif médical

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Date de fabrication

Reference and designation of the articles concerned:

MBR Mini Adjustable Connecting Rod Herbst Type (Pair)

Female part (upper) Male (lower) part

You have purchased our Mini-connecting rod system and we thank you for it. **Read this notice carefully.** The materials used for manufacturing have a perfect biological tolerance validated within the framework of the regulations relating to medical devices. Sold clean but not sterile, it is recommended, before use, to apply the cleaning rules in force in laboratories and dental offices. Please observe the storage instructions.

Intended use: intended for the orthodontic treatment of malocclusions in children and adults

Technical indication:

Used to make an orthodontic appliance consisting of two lap aligners. The mini-connecting rod is fixed on the gutter on plates inserted in the resin.

Clinical indication:

Orthodontic treatment of Class II malocclusions in growing children.

Clinical benefit:

Obtaining Class I

Warning:

Use is limited in orthodontic treatment. Use in patients is the responsibility of the practitioner.

Contraindication:

Not suitable for people with allergies to Nickel and Cobalt

Laboratory

ASSEMBLY: It is recommended that you have read the manufacturing protocol with PUL Concept

The assembly consists of connecting the lower gutter to the upper aligner (previously manufactured according to the PUL Concept protocol) with the Mini-connecting rod system (ref MBR) and platinum fasteners (ref APV) on the articulator in the corrected position (according to the practitioner's specifications).

Put the screw in the part of the system provided for this purpose and insert the assembly into the plate in the gutter; then tighten the screw using the mounting wrench (ref TCH1.5).

Practitioners

ACTIVATION: It is recommended to have attended clinical training on the PUL concept with Dr. CALLABE

The activation of the Mini-connecting rod system is done by unscrewing the nut mounted on the thread of the lower male part using the activation key (ref CAMB); 3 turns correspond to 1mm of activation.

Patients

DEVICE PORT:

The device should be worn for 12 to 14 hours a day. The device should not be used for eating, brushing teeth, playing wind music (flute, trumpet, etc.), and for sports where there is a risk of breakage or loss.

For the activation of the retractors, follow the instructions of the practitioner

CLEANING:

Before and after each use, clean the appliance under water with a toothbrush (different from the one for teeth) to avoid tartar.



PRECAUTIONS:

When removing the device, it should always be placed in its box to protect it. Never put it in a tissue or paper towel, it may be thrown in the trash. In case of cheek injury, slide a rolled compress between the cheek and the metal system. Do not put the device within reach of a dog, he may bite and break it.

Any serious incident that has occurred in connection with a medical device in the Pul Concept® range must be notified immediately to your healthcare professional and to PUL Concept on +33 2 54 75 54 80 Email: pulconcept@pulconcept.com; as well as to the competent health authority in your country.

PUL Concept
ZA The Plants
7 rue Pierre et Marie Curie
41140 NOYERS SUR CHER – France
02.54.75.54.80
pulconcept@pulconcept.com

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