

# DOCTOR \_\_\_\_\_ PRESCRIPTION PUL

Name:

Age:

Return date:

File number:

**TYPE OF PUL :**      PUL1      PUL2      PUL3      PULW      PUL Twin      PUL Snorlight (anti-snoring)

## **UPPER:**

- Simple (without screw, nor loop) - expansion screw - 3D expansion screw
- TMA Coffin loop - stainless steel Coffin loop
- Frame: mesially to 16/26 - mesially to 14/24
- Anterior hooks for elastic - anterior metal arch - anterior resin band – without anterior band
- Ball claps hooks: between 4&5 - between 5&6
- Schwartz spring (to push forward the incisors): TMA - stainless steel
- Other indications: \_\_\_\_\_

## **LOWER:**

- Simple (without screw, nor loop) - expansion screw - TMA lingual arch with 2 omegas
- Anterior metal arch - anterior resin band – without anterior band
- Grind lingual resin at 2 or 3 mm distance from the incisors (to upright them)
- Ball claps hooks: between 4&5 – between 5&6
- Other indications: \_\_\_\_\_

## **OCCLUSION WITH PUL:**

- Class I canines or premolars (in case of Class II division 1 or subdivision)
- Edge to edge (in case of Class II division 2 or Class III (PUL3 need bite wax))

## **PUL SNORLIGHT:**

- Mandibular advancement: 5 mm
- Recover occlusal resin on 33 to 43 (anti-bruxism, anti-clenching)
- Do not recover: free occlusal resin
- With anterior resin band